



Pour-thru cell kit for DR3900 and DR1900

Αρ. προϊόντος:

LZV899

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Permits rapid colorimetric measurements without requiring individual glass sample cells. The Pour-Thru Cell fits into the cell compartment, and solutions being analyzed are poured into the cell by means of a funnel and tubing. After measurement, solution from the cell flows through tubing to a drain. The Pour-Thru Cell Kit provides improved measurement accuracy and convenience.