



## Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Αρ. προϊόντος:

2745050

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Chlorine exists in water as "free" chlorine (hypochlorite) or "combined" chlorine (chloramines). These chlorine test strips streamline water testing for field technicians, enabling swift analysis of multiple samples within minutes and facilitating on-site evaluations without the need for extensive preparation or cleanup. Widely utilized in laboratories globally, our chlorine test strips serve as invaluable tools for pre-test screening, swiftly detecting potential interferences with lab testing procedures. Each package contains 50 tests.

---

### Προδιαγραφές

Εύρος μέτρησης :	0 - 10 mg/L $\text{Cl}_2$
Κατηγορία:	Test Strip
Μικρότερα βήματα ανάλυσης :	Steps: 0, 0.5, 1, 2, 4, 10
Μονάδες μέτρησης:	50 tests
Μοντέλο σετ μέτρησης :	TS
Παράμετρος :	Free and Total Chlorine
Υποσημείωση :	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$ ; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.