



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Αρ. προϊόντος:

2890200

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Προδιαγραφές

Εύρος μέτρησης :	0 - 600 mg/L Cl_2
Κατηγορία:	Test Strip
Μικρότερα βήματα ανάλυσης :	Steps: 0, 25, 50, 100, 200, 400, 600
Μονάδες μέτρησης:	100 tests
Παράμετρος :	Free Chlorine, High Range
Υποσημείωση :	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.