



Bottle, sampling, presterilized, round, 50/pk

Αρ. προϊόντος: 2495050

EUR Τιμή: Επικοινωνήστε μαζί μας
Δεν εμφανίζεται η ημερομηνία αποστολής

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen Medium.

Προδιαγραφές

Dechlorinating Agent?:	No
Όγκος/Μέγεθος πακέτου :	50 Stück