



## Sension+ pH gel electrode with ceramic diaphragm (general use)

Αρ. προϊόντος:

LZW5050T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Low cost gel-filled pH electrode with ceramic diaphragm. Automatic temperature compensation. Polycarbonate body. Range 0-14 pH, 0-80 °C. 1 m cable. For portable meters.

This low maintenance pH electrode is suitable for measurements in aqueous media in general. A storage protector prevents breakage and avoids over pressurisation of the diaphragm.

**Ergonomic and robust**

**Integrated temperature sensor**

**Low maintenance**

---

### Προδιαγραφές

|                      |                            |
|----------------------|----------------------------|
| Accuracy:            | 0.02 pH with pH1 Meter     |
| Electrode junction : | Ceramic pin                |
| Διάλυμα πλήρωσης :   | Non-refillable gel         |
| Διαστάσεις (BxM):    | 12 mm x 85 mm              |
| Εγγύηση :            | 6 μήνες                    |
| Εύρος θερμοκρασίας : | 0 - 80 °C                  |
| Εύρος μέτρησης :     | 0 - 14 pH                  |
| Θερμική αντίσταση :  | Pt1000                     |
| Μήκος :              | 85 mm                      |
| Μήκος καλωδίου :     | 1 m                        |
| Παράμετρος :         | pH                         |
| Σύνδεσμος :          | Sension+ MP5               |
| Τύπος :              | Φορητό                     |
| Τύπος αισθητηρίου :  | Glass                      |
| Τύπος ηλεκτροδίων :  | Portable                   |
| Υλικό :              | Sensor Body: Polycarbonate |