



Be Right™



## Hardness, iron, and pH test kit, model HA-62

Αρ. προϊόντος:

183700

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Model HA-62

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

**Popular tests at a great price in Hach's convenient and protective test kit carrying case.**

**Uses color disc and drop count titration methods**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Προδιαγραφές

Αριθμός δοκιμών :

100

Εύρος μέτρησης :

6.6 - 8.4 pH units

Εφαρμογή:

Water Quality

Κατηγορία:

Χρωματικός δίσκος

Μοντέλο :

HA-62

Παράμετρος :

Hardness, Iron, and pH