



Be Right™



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Αρ. προϊόντος:

223002

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Προδιαγραφές

Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 3,4 mg/L
Εφαρμογή:	Water Quality
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	CN-39WR
Παράμετρος :	Chlorine, Hardness, Iron, and pH