



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Αρ. προϊόντος:

143601

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Προδιαγραφές

Smallest increment:	1.25
	2
	5
Αριθμός δοκιμών :	200
Εύρος μέτρησης :	5 - 100 mg/L CO ₂
Εύρος μέτρησης (2):	1.25 - 25 mg/L
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	CA-23
Όνομα μεθόδου :	Drop count titration/ Sodium hydroxide Phenolphthalein Carbon Dioxide - As CO ₂
Παράμετρος :	D
Τύπος θήκης :	
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.