



Test kit, total chlorine, model CN-DT, 20-2000 mg/L

Αρ. προϊόντος:

2471100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Kit contains all required reagents and apparatus to perform analysis in the field

Προδιαγραφές

Smallest increment:	0.2 - 5.0
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	20 - 2,000 mg/L Cl ₂
Κατηγορία:	Ψηφιακός τιτλοδότης
Μοντέλο σετ μέτρησης :	CN-DT
Όνομα μεθόδου :	Digital Titrator/Thiosulfate
Παράμετρος :	Chlorine, mid and high range - As Total Cl ₂
Τύπος θήκης :	F
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.