



Dissolved oxygen test kit, model OX-2P

Αρ. προϊόντος:

146900

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Προδιαγραφές

Smallest increment:	0.2
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	OX-2P
Όνομα μεθόδου :	Drop count titration/ Modified Winkler
Παράμετρος :	Dissolved Oxygen - As DO
Τύπος θήκης :	D