



Test kit, alkalinity, model AL-DT

Αρ. προϊόντος:

2063700

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Αριθμός δοκιμών :	100
Εύρος μέτρησης :	10 - 4,000 mg/L CaCO ₃
Κατηγορία:	Ψηφιακός τιτλοδότης
Μικρότερα βήματα ανάλυσης :	0.1-10
Μοντέλο σετ μέτρησης :	AL-DT
Όνομα μεθόδου :	Digital Titrator/Sulfuric acid
Παράμετρος :	Alkalinity - (P) & Total (MO)
	as CaCO ₃
Τύπος θήκης :	F
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.