



Be Right™



Test kit, alkalinity, model AL-AP

Αρ. προϊόντος:

2444300

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:

0.4 gpg

1 gpg

Αριθμός δοκιμών :

100

Εύρος μέτρησης :

0 - 8 gpg & 1 - 20 gpg

Κατηγορία:

Drop Count Titration

Μοντέλο σετ μέτρησης :

AL-AP

Όνομα μεθόδου :

Drop count titration/Sulfuric acid

Παράμετρος :

Alkalinity - (P) & Total (MO) As CaCO_3

Τύπος θήκης :

B

Υποσημείωση :

*mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange