



Alkalinity test kit AL-AP, 5 - 100 / 20 - 400 mg/L CaCO₃ 100 pcs

Αρ. προϊόντος:

2444301

Επικίνδυνο

EUR Τιμή:
Διαθέσιμο

Επικοινωνήστε μαζί μας



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Προδιαγραφές

Smallest increment:	5
	20
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	5 - 100 mg/L & 20 - 400 mg/L
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	AL-AP, MG-L
Όνομα μεθόδου :	Drop count titration/Sulfuric acid
Παράμετρος :	Alkalinity - (P) & Total (MO) As CaCO ₃
Τύπος θήκης :	B
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.