



Sension+ 5014T pH liquid combination electrode, with silver ion barrier

Αρ. προϊόντος:

LZW5014T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

pH liquid electrode (refillable) with 2 ceramic diaphragms and silver ion barrier. Automatic temperature compensation. Glass body. Range 0-14 pH, -10 to 100 °C. 1 m cable. For benchtop meter.

This pH electrode is suitable for aqueous media in general and for samples containing sulphurs, reducing sugars or other substances which react with silver ions.

Reference electrode stability due to sleeved cartridge and Ag⁺ ion barrier design.

For samples containing sulphurs, reducing sugars or other substances which react with silver ions.

Ag⁺ ion barrier reference design.

High stability and long lifetime

Encapsulated reference system.

Προδιαγραφές

Accuracy:	0.002 pH with pH31 Meter
Electrode junction :	2 x Ceramic pin
Διάλυμα πλήρωσης :	LZW9500.99 KCl 3M
Διαστάσεις (BxM):	12 mm x 130 mm
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	TRIS compatible
Εύρος θερμοκρασίας :	-10 - 100 °C
Εύρος μέτρησης :	0 - 14 pH
Θερμική αντίσταση :	Pt1000
Μήκος :	130 mm
Μήκος καλωδίου :	1 m
Παράμετρος :	pH
Σύνδεσμος :	BNC for pH + Banana for Pt 1000
Τύπος αισθητηρίου :	Glass
Τύπος ηλεκτροδίων :	Laboratory
Υλικό :	Sensor Body: Glass