



Be Right™



## Sension+ pH combination electrode 5021T with A.T.C. for difficult samples

Αρ. προϊόντος:

LZW5021T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

pH combination electrode 5021T with sleeve diaphragm, refillable electrolyte and automatic temperature compensation. Glass body. Range 0-14 pH, 0-60 °C. 1 m cable. For benchtop instruments.

This electrode is suitable for the most difficult samples, e.g. for low ionic strength samples, creams, emulsions, very viscous samples, wines, paint and samples with colloids. The encapsulated reference ensures high stability and long lifetime.

### Fast responding

Large cylindrical pH glass membrane, sleeve diaphragm.

### High stability and long lifetime

Encapsulated Ag/AgCl reference system

### Automatic temperature compensation

Built-in Pt1000 temperature sensor.

## Προδιαγραφές

Accuracy:	0.002 pH with pH31 Meter
Electrode junction :	Open with sleeve
Διάλυμα πλήρωσης :	LZW9500.99
Διαστάσεις (BxM):	12 mm x 130 mm
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	The encapsulated reference and the clog-free sleeve junction ensure high stability and long lifetime. The large cylindrical pH membrane provides fast response.
Εύρος θερμοκρασίας :	0 - 60 °C
Εύρος μέτρησης :	0 - 14 pH
Θερμική αντίσταση :	Pt1000
Μήκος :	130 mm
Μήκος καλωδίου :	1 m
Παράμετρος :	pH
Σύνδεσμος :	BNC for pH + Banana for Pt 1000
Τύπος αισθητηρίου :	Glass
Τύπος ηλεκτροδίων :	Laboratory

Υλικό :

Sensor Body: Glass