



Sension+ Multisensor for pH, ORP & conductivity

Αρ. προϊόντος:

LZW5048.97.0002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Multisensor for pH (0-14), ORP (± 2000 mV), conductivity (25-30000 $\mu\text{S}/\text{cm}$) & temperature (0-80 °C). Automatic temperature compensation. Cable length 1.5m. Specific for Sension+ MM150 multimeter.

With Polycarbonate body, gel electrolyte and ceramic diaphragm.

Attachable calibration tubes can screw into handle of electrode (only 2 mL of buffer or standard). Storage protector prevents breakage and avoids over pressurization of the diaphragm.

Accurate

Easy to use

Rugged & durable

Προδιαγραφές

Accuracy:	1 mV with MM150 Meter
Electrode junction :	Ceramic
Διάλυμα πλήρωσης :	Non-refillable gel
Διαστάσεις (BxM):	12 mm x 85 mm
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	pH, Conductivity, ORP
Εύρος θερμοκρασίας :	0 - 80 °C
Εύρος μέτρησης :	Conductivity: 20 $\mu\text{S}/\text{cm}$ - 200 mS/cm
Εύρος μέτρησης (2):	pH: 0 - 14 pH
Θερμική αντίσταση :	Pt1000
Μήκος :	85 mm
Μήκος καλωδίου :	1.5 m
Παράμετρος :	pH, ORP, Conductivity
Σύνδεσμος :	Sension+ MP8
Τύπος :	Φορητό
Τύπος αισθητηρίου :	pH: Glass
	ORP: Platinum

Τύπος ηλεκτροδίων :

Conductivity: Platinum

Υλικό :

Portable combination pH, ORP, conductivity

Sensor Body: Polycarbonate