



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Αρ. προϊόντος:

108076

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Προδιαγραφές

Διάμετρος, εξωτερική :	83.6 mm with spout
Μονάδες μέτρησης:	6 /pk
Σώμα υλικού:	Polypropylene
Υλικό :	Polypropylene, Low Form
Χωρητικότητα :	250 mL