



Filter, Membrane, Pore Size 0.45 µm, Diameter 47 mm, 1000/ pk, without pad

Αρ. προϊόντος:

1353002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

These 47 mm diameter, 0.45 µm pore size GN-Metricel membrane filters meet or surpass all U.S. Environmental Protection Agency (USEPA) and American Society for Testing and Materials (ASTM) specifications to ensure complete retention and maximum recovery of bacteria. Black membranes for yeast and mold are also available. Black membranes make yeast and mold colonies more distinct and easier to count. Filters are available in convenient, presterilized, individually sealed packs (P/N 13530-01 and 24903-00) or in resealable envelopes with 10 filters and 10 absorbent pads each (P/N 21269-00). Phosphate free.

Ensure complete retention and maximum recovery of bacteria

Phosphate free

Προδιαγραφές

Filter Type:	Filter
Διάμετρος :	47 mm
Μέγεθος πόρου :	0.45 µm
Όγκος/Μέγεθος πακέτου :	1000 /pk
Υλικό :	Mixed-Cellulose Ester