



Cyanuric acid test kit, model CY-3

Αρ. προϊόντος:

185102

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:	20
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	20 - 100 mg/L
Εφαρμογή:	Pool and Spa
Μοντέλο :	CY-3
Μοντέλο σετ μέτρησης :	CY-3
Όνομα μεθόδου :	Turbidimetry
Παράμετρος :	Cyanuric Acid
Τύπος θήκης :	D