



Test kit, hardness & iron color disc, model HA-77

Αρ. προϊόντος:

202300

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

The kit to use for more resolution in a visual test

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Εφαρμογή:	Water Quality
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	HA-77
Παράμετρος :	Hardness and Iron