



Test kit, chloride, model CD-51

Αρ. προϊόντος:

208601

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:	500
	5000
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	500 - 10000 mg/L & 5000 - 100000 mg/L Cl ⁻
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	CD-51
Όνομα μεθόδου :	Drop count titration/
	Silver nitrate
Παράμετρος :	Chloride - As Cl ⁻
Τύπος θήκης :	B
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.