



Test kit, formaldehyde, model FM-1

Αρ. προϊόντος:

2183100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:	0.05%
	0.5%
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0,05 - 1% & 0,5 - 10%
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	FM-1
Όνομα μεθόδου :	Drop count titration/ Thymolphthalein
Παράμετρος :	Formaldehyde - As CH_2O
Τύπος θήκης :	D