



## Test kit, chromium, model CH-14

Αρ. προϊόντος:

222702

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Προδιαγραφές

Smallest increment:

5

50

Αριθμός δοκιμών :

100

Εύρος μέτρησης :

5 - 100 mg/L & 50 - 1000 mg/L Cr<sup>6+</sup>

Κατηγορία:

Drop Count Titration

Μοντέλο σετ μέτρησης :

CH-14

Όνομα μεθόδου :

Drop count titration/

Thiosulfate

Παράμετρος :

Chromium, high range - As Hexavalent Cr

Τύπος θήκης :

D