



Test kit, chlorine (total), model CN-65

Αρ. προϊόντος:

225401

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:	0.2
	1
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0,2 - 4,0 mg/L & 1 - 20 mg/L
Κατηγορία:	Drop Count Titration
Μικρότερα βήματα ανάλυσης :	1
Μοντέλο σετ μέτρησης :	CN-65
Όνομα μεθόδου :	Drop count titration/ Thiosulfate
Παράμετρος :	Chlorine, mid and high range - As Total Cl ₂
Τύπος θήκης :	D