



Test kit, ascorbic acid, model ASC-1

Αρ. προϊόντος:

2308100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Smallest increment:	10
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	10 - 200 mg/L
Κατηγορία:	Drop Count Titration
Μοντέλο σετ μέτρησης :	ASC-1
Όνομα μεθόδου :	Drop count titration/
	Iodometric
Παράμετρος :	Ascorbic Acid
Τύπος θήκης :	D
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.