



Microfil funnel & filter, 100 mL, 0.45 µm 47 mm white gridded (millipore)

Αρ. προϊόντος:

2586300

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 150 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45 µm pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

Προδιαγραφές

| | |
|-------------------------|-----------------------|
| Διάμετρος : | 47 mm |
| Μέγεθος πόρου : | 0.45 µm |
| Όγκος/Μέγεθος πακέτου : | 150 |
| Υλικό : | Mixed-Cellulose Ester |
| Χωρητικότητα : | 100 mL |