



Test kit color disc chlorine dioxide, model CLD-2

Αρ. προϊόντος:

2842800

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Monitor ClO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Προδιαγραφές

Parameter:	Chlorine Dioxide
Smallest increment:	0.1
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 2.0 ppm ClO ₂
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	CLD-2
Μοντέλο σετ μέτρησης :	CLD-2 ²
Όνομα μεθόδου :	Color Disc/DPD
Παράμετρος :	Chlorine Dioxide
Τυπικό εύρος μέτρησης:	0 - 2.0 ppm ClO ₂
Τύπος θήκης :	D