



Brush, buret, 50 mL

Αρ. προϊόντος:

68600

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Ideally shaped for quick and thorough cleaning of burets and other similarly shaped vessels. Unbleached white nylon bristles are supported by twisted galvanized steel wire. Brush portion measures 8 cm (3.1") in length. Fan style tip. 2.0cm (0.79") Dia. 91cm (35.8") L

Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

1 each