



## Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Αρ. προϊόντος:

2966700

### ΚΑΤΑΡΓΗΜΕΝΟ ΣΤΟΙΧΕΙΟ

Αυτό το στοιχείο δεν είναι πλέον διαθέσιμο.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

---

### Προδιαγραφές

Μέγεθος πόρου :	10 µm
Όγκος/Μέγεθος πακέτου :	25 /pk
Υλικό :	Cellulose