



Maintenance Kit, CL17 Chlorine Analyzer (for instruments sold before April 2000)

Αρ. προϊόντος:

4721600

ΚΑΤΑΡΓΗΜΕΝΟ ΣΤΟΙΧΕΙΟ

Αυτό το στοιχείο δεν είναι πλέον διαθέσιμο.

CL17 maintenance kit, for instruments purchased before April 2000. Kit includes reagent tubing, colorimeter lamp, reagent caps and fittings for annual replacement (pump module tubing to be replaced at 3-6 month interval). For instruments purchased after March 2000, order PN 5444301.