



Be Right™



C501X Cable SubD-9f/2m/SubD-9f

Αρ. προϊόντος:

A95X501

EUR Τιμή:
Διαθέσιμο

Επικοινωνήστε μαζί μας

RS232 Serial Connection Cable, SubD-9

The CCTO is a pyrex glass body cell designed for testing corrosion on large flat samples. The shape of the cell combined to the large area of the working electrode and the stainless steel counter electrode make it particularly suitable for testing low conductivity media and/or highly resistive samples. The setup includes a calomel reference electrode mounted on a vertically adjustable Luggin capillary designed to minimize ohmic drop.