



Test Kit, chlorine (free & total) & pH, model CN-67

Αρ. προϊόντος:

1411100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Αριθμός δοκιμών :	300
Εύρος μέτρησης :	6.6 - 8.4 pH units
Εύρος μέτρησης (2):	0 - 3,4 mg/L
Εφαρμογή:	Pool and Spa
Κατηγορία:	Χρωματικός δίσκος