



Deluxe water conditioning demonstration test kit, model AD-16E

Αρ. προϊόντος: 148410

EUR Τιμή: Επικοινωνήστε μαζί μας
Τηλεφωνήστε για την ημερομηνία αποστολής

Επικίνδυνο



Designed to boost your testing capabilities, Hach's most complete water conditioning test kit includes tests and demonstrations for seven critical water quality parameters. A optional softener regeneration test kit is available (P/N - 2666000).

Uses color disc and drop count titration methods

In addition to soap consumption and softener demonstrations, includes tests to measure iron, pH, and nitrate

Also includes the Total Dissolved Solids (TDS) Pocket Pal™ Tester

Comes complete with a sturdy carrying case and features viewing tube, water softener, and a compartmentalized reagent cabinet to keep chemistries neat and well-organized

Προδιαγραφές

Εφαρμογή:	Demonstration
Κατηγορία:	Χρωματικός δίσκος