



Testkit, glycol visual, model EG-1

Αρ. προϊόντος:

2186400

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Αριθμός δοκιμών :	25
Εύρος μέτρησης :	Presence or absence of glycol, 0- greater than 150 mg/L
Μοντέλο σετ μέτρησης :	EG-1
Όνομα μεθόδου :	Visual/ Appearance of color
Παράμετρος :	Glycol - Ethylene glycol in oil or water
Τύπος θήκης :	G