



Be Right™



Test kit, manganese color disc, model MN-PAN

Αρ. προϊόντος:

2350800

Επικίνδυνο

EUR Tμη:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

Προδιαγραφές

Parameter:	Manganese
Smallest increment:	0.05
Αριθμός δοκιμών :	50
Εύρος μέτρησης :	0 - 0,70 mg/L Mn
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	MN-PAN
Μοντέλο σετ μέτρησης :	MN-PAN
Όνομα μεθόδου :	Color disc/PAN
Παράμετρος :	Manganese - As Mn
Τυπικό εύρος μέτρησης:	0 - 0.70 mg/L Mn
Τύπος θήκης :	D