



Test kit cyanide, model CYN-3

Αρ. προϊόντος:

201002

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Προδιαγραφές

Parameter:	Cyanide
Smallest increment:	0.02
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 0,3 mg/L CN ⁻
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	CYN-3
Μοντέλο σετ μέτρησης :	CYN-3
Όνομα μεθόδου :	Color disc/ Pyridine-Pyrazolone
Παράμετρος :	Cyanide - As free CN ⁻
Τυπικό εύρος μέτρησης:	0 - 0.3 mg/L CN ⁻
Τύπος θήκης :	R