



Test kit, sulfate, model SF-1

Αρ. προϊόντος:

225100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Προδιαγραφές

Smallest increment:	50
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	50 - 200 mg/L
Μοντέλο σετ μέτρησης :	SF-1
Όνομα μεθόδου :	Extinction/Turbidimetric
Παράμετρος :	Sulfate - As SO_4^{2-}
Τύπος θήκης :	D