



Be Right™



## Test kit, alkalinity, model AL-AP

Αρ. προϊόντος: 2444300  
EUR Τιμή: Επικοινωνήστε μαζί μας  
Αποστολή εντός 2 εβδομάδων

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Προδιαγραφές

Smallest increment:	0.4 gpg
	1 gpg
Αριθμός δοκιμών :	100
Βάρος μεταφοράς :	1
Εύρος μέτρησης :	- 8 gpg
Εύρος μέτρησης (2):	1 - 20 gpg**
Κατηγορία:	Ογκομετρική τιτλοδότηση
Μοντέλο σετ μέτρησης :	AL-AP
Όνομα μεθόδου :	Drop count titration/Sulfuric acid
Παράμετρος :	Alkalinity - (P) & Total (MO) as CaCO <sub>3</sub>
Τύπος θήκης :	B
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.

\*\*1 gpg = 17.1 mg/L (MO) = Methyl Orange