



## Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

Αρ. προϊόντος:

224000

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as  $\text{NO}_2^-$  and contains 100 tests. Range (mg/L): 0-8/0-1600.

**Simple, low-cost colour disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Προδιαγραφές

Parameter:	Nitrogen, Nitrite
Smallest increment:	2
	40
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	Single Range
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο σετ μέτρησης :	NI-6
Όνομα μεθόδου :	Colour disc/
	Ferrous sulphate
Παράμετρος :	Nitrogen, Nitrite, high range - as $\text{NO}_2^-$
Τι περιέχει το καλάθι:	Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.
Τύπος θήκης :	D

### Τι περιέχει το καλάθι

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.