



Beaker, polypropylene, low form, 50 mL, 12/pk

Αρ. προϊόντος:

108071

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Προδιαγραφές

Διάμετρος, εξωτερική :

51.1 mm with spout

Μονάδες μέτρησης:

12 /pk

Σώμα υλικού:

Polypropylene

Υλικό :

Polypropylene, Low Form

Χωρητικότητα :

50 mL