



Λογισμικό εφαρμογής για ενζυμική ανάλυση τροφίμων

Αρ. προϊόντος:

LZW943

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

The application software LZW937 is a compilation of enzymology tests

manufactured by R-Biopharm AG, Darmstadt that can be performed with the spectrophotometer DR6000 and carousel insert LZW902.99.00001. The methods can be modified for individual settings.

Methods are ready to go.

The methods are convenient to upload via USB.

Keep control of time depending multiple measurements.

A timer function alerts about the next step and starts evaluation of multiple samples automatically.

Get the final result immediately after the analysis.

The photometer does the complex result calculation for you.

Απαιτούμενα παρελκόμενα

- Περιστρεφόμενη βάση 7 x 1 cm (Item LZW902.99.00002)