



Sension+ pH combination electrode (dirty samples)

Αρ. προϊόντος:

LZW5011T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

pH electrode with open circular diaphragm and solid electrolyte. Automatic temperature compensation. Glass body. Range 2-14 pH, 0-80 °C. 1 m cable. For benchtop meter.

The electrode is suitable for waste water, dirty samples, viscous media, and food applications. The open diaphragm around the membrane facilitates the measuring in some media where other electrodes have certain difficulties.

Automatic temperature compensation

Built-in Pt1000 temperature sensor.

Προδιαγραφές

Accuracy:	0.002 pH with pH31 Meter
Electrode junction :	Open
Διάλυμα πλήρωσης :	Non-refillable solid polymer
Διαστάσεις (BxM):	12 mm x 130 mm
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	Solid gel for high solids content
Εύρος θερμοκρασίας :	0 - 80 °C
Εύρος μέτρησης :	2 - 14 pH
Θερμική αντίσταση :	Pt1000
Μήκος :	130 mm
Μήκος καλωδίου :	1 m
Παράμετρος :	pH
Σύνδεσμος :	BNC for pH + Banana for Pt 1000
Τύπος αισθητηρίου :	Glass
Τύπος ηλεκτροδίων :	Laboratory
Υλικό :	Sensor Body: Glass