



Sension+ 5130 Portable Polarographic Dissolved Oxygen Probe

Αρ. προϊόντος:

LZW5130.97.0002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

DO polarographic probe with 1 meter cable and MP5 connector for use with portable meters.

Polarographic dissolved oxygen (DO) probe with a durable plastic body, replaceable oxygen-permeable membrane, refillable reference electrolyte, and built-in temperature sensor. The probe has a fixed 1 meter cable and MP5 connector dedicated for use with Portable Dissolved Oxygen meters. Probe is shipped with 30 mL bottle of electrolyte solution, two membrane modules and polishing strips.

The electrode is ideal for dissolved oxygen measurement in general aqueous applications – particularly in a field environment.

Note: Sension+ 5130 DO Probe is compatible with versions of the Sension+ DO6 meter produced since January 2013. To ensure compatibility with your Sension+ DO6 meter, please contact Support for assistance.

Προδιαγραφές

Accuracy:	±0.5 mg/L
Διάλυμα πλήρωσης :	2759123
Διαστάσεις (BxM):	12 mm x 120 mm (χωρίς χειρολαβή)
Εγγύηση :	6 μήνες
Εύρος θερμοκρασίας :	0 - 50 °C
Εύρος μέτρησης :	0.03 mg/L - saturation
Θερμική αντίσταση :	Integrated 30kohm NTC
Μήκος :	120 mm (χωρίς χειρολαβή)
Μήκος καλωδίου :	1 m
Παράμετρος :	Dissolved Oxygen
Σύνδεσμος :	Sension+ MP5
Τύπος :	Φορητό
Τύπος αισθητήριου :	Replaceable PTFE membrane; Pt cathode; Ag anode
Τύπος ηλεκτροδίων :	Portable Polarographic Dissolved Oxygen
Υλικό :	Probe body: ABS & POM (nylon)

Απαιτούμενα παρελκόμενα

- Διάλυμα πλήρωσης, εσωτερικό, πολαρογραφικό DO, 50 mL (Item 2759126)