



Sension+ pH combination electrode with open diaphragm (dirty samples)

Αρ. προϊόντος:

LZW5051T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

pH electrode with open circular diaphragm and solid electrolyte. Automatic temperature compensation. Glass body. Range 2-14 pH, 0-80 °C. 1m cable. For portable meters.

The electrode is suitable for (very dirty) waste water, milk, juices, and food applications. A polypropylene working protector prevents accidental breakage of the electrode.

Automatic temperature compensation

Built-in Pt1000 temperature sensor.

Ergonomic and robust

Low maintenance

Προδιαγραφές

Accuracy:	0.02 pH with pH1 Meter
Electrode junction :	Open
Διάλυμα πλήρωσης :	Non-refillable solid polymer
Διαστάσεις (BxM):	12 mm x 85 mm (χωρίς χειρολαβή)
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	With Working Protector. Solid gel for high solids content.
Εύρος θερμοκρασίας :	0 - 80 °C
Εύρος μέτρησης :	2 - 14 pH
Θερμική αντίσταση :	Pt1000
Μήκος :	85 mm (χωρίς χειρολαβή)
Μήκος καλωδίου :	1 m
Παράμετρος :	pH
Σύνδεσμος :	Sension+ MP5
Τύπος :	Φορητό
Τύπος αισθητηρίου :	Glass
Τύπος ηλεκτροδίων :	Portable
Υλικό :	Sensor Body: Glass