



## Sension+ pH combination electrode (difficult samples)

Αρ. προϊόντος:

LZW5052T.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Gel-filled pH combination electrode with annular porous PTFE diaphragm and automatic temperature compensation Pt 1000. Glass body. Range 0-14 pH, 0-100 °C. 1m cable, MP5 connector, for portable meters.

The electrode is suitable for difficult media, e.g. for low ionic strength samples, creams or viscous samples. Measurements at high temperature up to 100 °C. Supplied with a working protector to avoid accidental breakage.

### Automatic temperature compensation

Built-in Pt1000 temperature sensor.

### Ergonomic and robust

### Low maintenance

---

## Προδιαγραφές

|                         |                                                                         |
|-------------------------|-------------------------------------------------------------------------|
| Accuracy:               | 0.02 pH with pH1 Meter                                                  |
| Electrode junction :    | Porous annular PTFE                                                     |
| Διάλυμα πλήρωσης :      | Non-refillable gel                                                      |
| Διαστάσεις (BxM):       | 12 mm x 85 mm (χωρίς χειρολαβή)                                         |
| Εγγύηση :               | 6 μήνες                                                                 |
| Ειδικά χαρακτηριστικά : | With Working Protector. For low ionic strength and/or high temperature. |
| Εύρος θερμοκρασίας :    | 0 - 100 °C                                                              |
| Εύρος μέτρησης :        | 0 - 14 pH                                                               |
| Θερμική αντίσταση :     | Pt1000                                                                  |
| Μήκος :                 | 85 mm (χωρίς χειρολαβή)                                                 |
| Μήκος καλωδίου :        | 1 m                                                                     |
| Παράμετρος :            | pH                                                                      |
| Σύνδεσμος :             | Sension+ MP5                                                            |
| Τύπος :                 | Φορητό                                                                  |
| Τύπος αισθητηρίου :     | Glass                                                                   |
| Τύπος ηλεκτροδίων :     | Portable                                                                |
| Υλικό :                 | Sensor Body: Glass                                                      |