



PHC2005 Robust combined pH electrode, Red Rod, BNC

Αρ. προϊόντος:

E16M500

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: porous pin. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Προδιαγραφές

Εγγύηση : 12 μήνες

Τύπος ηλεκτροδίου : PHC2005