



PHC2085 Combined pH electrode, Red Rod, Temp. sensor

Αρ. προϊόντος:

E16M501

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: porous pin. With temperature sensor.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Προδιαγραφές

Electrode junction :	Porous pin
Διάλυμα πλήρωσης :	Refillable: Sat. KCl
Διαστάσεις (BxM):	12 mm x 103 mm
Εγγύηση :	12 μήνες
Ειδικά χαρακτηριστικά :	Red-Rod reference system
Εύρος θερμοκρασίας :	-10 - 100 °C
Εύρος μέτρησης :	0 - 12 pH
Θερμική αντίσταση :	NTC 10 kohm
Μήκος :	103 mm
Μήκος καλωδίου :	1 m
Τύπος αισθητηρίου :	Glass
Τύπος ηλεκτροδίων :	Laboratory
Υλικό :	Sensor body: Epoxy