



Tensette Pipette 1.0 - 10.0 mL

Αρ. προϊόντος:

1970010

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Perform routine laboratory pipetting quickly, economically, and accurately with Hach's automatic Tensette Pipette. Two models dispense from 0.1 - 1 mL or from 1 - 10 mL in 0.1 or 1.0 mL increments, respectively. The pipette saves time, reduces glassware cleanup and assures accuracy. A Set includes one pipette, a manual and 50 disposable tips.

Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

1 each