



REF201 Universal reference electrode, 7.5mm, Red Rod

Αρ. προϊόντος:

E21M009

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

General-purpose Red Rod reference electrode with banana plug. Liquid junction: porous pin, salt bridge solution: sat. KCl. Temperature range -10 to 100°C.

Red Rod electrodes feature:

- Reproducible measuring conditions
- Fast and accurate response
- Wide temperature range
- Stable signal
- Trouble-free liquid junction means less risk of clogging

Προδιαγραφές

Εγγύηση : 12 μήνες

ηλεκτροδίων αναφοράς : REF201