



## GK ANNEX Maintenance Kit for Electrodes, Complete (Radiometer Analytical)

Αρ. προϊόντος:

S91M001

### **ΚΑΤΑΡΓΗΜΕΝΟ ΣΤΟΙΧΕΙΟ**

Αυτό το στοιχείο δεν είναι πλέον διαθέσιμο.

For regular maintenance of combined pH, glass and reference electrodes with saturated potassium chloride salt bridge.

The GK ANNEX Maintenance Kit contains all you need for regular maintenance of combined pH, glass, and reference electrodes with saturated potassium chloride salt bridg. It contains: solutions for normal and intensive cleaning, KCl as saturated solution and crystals for refilling, GLP Logbook to keep note of calibration, and maintenance procedures, utensils and easy to follow instructions.

---

### Προδιαγραφές

Ειδικές σημειώσεις :

Kit contains: Solutions for normal and intensive cleaning, KCl as saturated solution and crystals for refilling, GLP logbook, Utensils and easy-to-follow instructions.

Παράμετρος :

GK ANNEX Electrode Maintenance Kit

Περιγραφή :

The Kit contains all you need for regular maintenance of combined pH, glass and reference electrodes with saturated potassium chloride salt bridge